

Table S2. Flow Cytometry Instrument Configuration

Antibody	Fluorochrome	Laser Emission Wavelength (nm)	Filter (ARIA IIIU)	Filter (SYNERGY)
CD3	Pacific-Blue	405	450/40	455/50
CD70	FITC	488	530/30	525/50
CD28	PE-CY5	488	695/40	665/30
CD28	PE-CY5	561	670/14	
CD40L	PE-Cy7	561	780/60	775/50
Kir	PE	561	582/15	577/25
CD11a	APC	633	660/20	665/30
CD11a	APC	R/O (594/685)		665/30
CD4	APC-Cy7	633	780/60	
CD4	APC-Cy7	R/O (594/685)		775/50

